

octo**flux**

## Procedures made easier with **SETS** designed for **CT INJECTORS**



### Patient safety is your priority...

- Spikes included in the CT sets avoid opening contrast medium vials: Asepsis is consequently improved without personnel injuries (cuts).
- Contamination risks are limited thanks to:
  - a drip chamber equipped with filters (air, fluid)
  - and non-return valves on the filling and injection set as well as patient line.

### Saving time and money...

- Limiting the number of operations:
  - syringe is filled without having to disconnect catheters or actuating stopcocks thanks to the non-return valves of the set,
  - the contrast medium is automatically administered to the patient during the injection phase.
- Contrast media is saved thanks to the low tube residual volume.
- Connections made easier: connections are made easier by the shape of the female Luer lock and spinlock male Luer lock.





## Filling and injection set:

Description	References	Packaging
<ul style="list-style-type: none"> <li>Length of tubing: 200 cm</li> <li>Maximum pressure: 300 PSI</li> <li>2 non-return valves (one on contrast line, one on patient line)</li> </ul> + a drip chamber for single-head injector	0170NF	25-unit box
<ul style="list-style-type: none"> <li>Length of tubings: patient line 40 cm and 10 cm contrast line 100 cm and 100 cm</li> <li>Maximum pressure: 300 PSI</li> <li>4 non-return valves (one on each contrast line, one on each patient line)</li> </ul> + 2 drip chambers for double-head injector	0170PS	25-unit box
<ul style="list-style-type: none"> <li>Length of tubings: patient line 150 cm and syringe to syringe 50 cm</li> <li>Maximum pressure: 300 PSI</li> </ul> Standard set, without non-return valve	0170NQ	25-unit box
<ul style="list-style-type: none"> <li>Length of tubings: patient line 150 cm and syringe to syringe 25 cm</li> <li>Maximum pressure: 300 PSI</li> </ul> 1 non-return valve on the saline line	0170RK	50-unit box
<ul style="list-style-type: none"> <li>Length of tubings: patient line 150 cm and syringe to syringe 25 cm</li> <li>Maximum pressure: 300 PSI</li> </ul> 2 non-return valves on the saline and the contrast lines	0170SR	50-unit box

## Patient line:

Description	References	Packaging
<ul style="list-style-type: none"> <li>Length of tubing: 20 cm</li> <li>Maximum pressure: 300 PSI</li> </ul> 1 non-return valve (on patient side)	0170NB	50-unit box
<ul style="list-style-type: none"> <li>Length of tubing: 40 cm</li> <li>Maximum pressure: 300 PSI</li> </ul> 1 non-return valve (on patient side)	0170NG	50-unit box
<ul style="list-style-type: none"> <li>Length of tubing: 60 cm</li> <li>Maximum pressure: 300 PSI</li> </ul> 1 non-return valve (on patient side)	0170NN	50-unit box
<ul style="list-style-type: none"> <li>Length of tubing: 120 cm</li> <li>Maximum pressure: 300 PSI</li> </ul> 1 non-return valve (on patient side)	0170NT	50-unit box